

# GEM Premier 4000

## Issues Log

<b>Hospital:</b>	<b>Ward:</b>
------------------	--------------

Instrument exterior and touchscreen should be cleaned and ready for use each day or as required.

### Cleaning Instructions

**Touch screen:** Clean the touchscreen using a soft cloth dampened slightly with water

**Exterior:** 1. Dampen a soft cleaning cloth with water or mild cleaning solution (make sure the cleaning cloth is only moist, not dripping wet).  
2. Remove any blood or dust for the outer surface of the instrument case using a clean, soft cloth moistened with 10% bleach.

Date	Time	Issue	Action Taken	Name